# Beer's Law Inquiry Activity

### Part 1: Colored Solutions and Spectrophotometers

#### **Introduction**

Take a minute to **explore** the simulation by clicking on the link provided. Specifically, check out the Beer's Law tab. (<u>http://phet.colorado.edu/sims/html/beers-law-lab/latest/beers-law-lab\_en.html</u>).

Using the virtual spectrophotometer in the Beer's Law tab, how does Concentration affect how much light is **absorbed** and **transmitted** through the solution?

Teacher initial:

#### **Investigating Absorption and Concentration**

- a. **Predict** what a graph of absorbance versus concentration would look like. Sketch your prediction below using a colored pencil.
- b. On the simulation, the solution chosen for you is of a generic "Drink Mix" and measure the Absorbance for 5 different concentrations, chosen for you in the table below, with the <u>fixed wavelength</u> setting and graph the data using a different colored pencil

Conc (M)	Conc (mM)	Abs
0		
0.05		
0.1		
0.2		
0.3		

				2		2	

Concentration

d. Based on Beer's Law (A = elC, A = absorbance, e = molar absorptivity, I = pathlength and C = concentration), would you expect the wavelength used to affect your absorbance versus concentration relationship?

Teacher initial: \_\_\_\_\_


## Investigating Absorption and Wavelength

3. Compare three solutions of different colors with the same pathlength.

		Preset Wavelength: Simulation default setting			Variable Wavelength: Set to same color as solution		
Solution name	Soln. Color	Beam Color	Value (nm)	Abs	Beam Color	Value (nm)	Abs

a. How are beam color, solution color, and absorbance related? What combinations give the most absorbance? Why?

b. Choose a solution and set concentration and pathlength. Keep these settings constant and graph the absorbance for 8 different wavelengths. Your absorbance range may vary depending on your solution settings.



d. Compare your absorbance spectrum sketch with a group that chose a different solution. Would you use the same wavelength of light to do spectroscopy experiments with different colored solutions? Why or why not?